



HALLOWELL EMC

000A6541

ICU Non-Rebreathing Module



Module shown with 2002Pro Model ventilator.

Converts your Pro Series Anesthesia Ventilator into an ICU Respirator and mounts conveniently under the vent. This unit was produced during COVID pandemic research and development and is a game changer for ICU care! A must-have for every facility that uses an HEMC Pro series Model Ventilator.



Panel Gauges show Air and O₂ psi from supply hoses.


This module easily attaches to a pedestal stand for better mobility.



Comes with

FiO₂ Approximation Guide & O₂ Blending Nomogram Card

- ◇ Equal flows yields about 60% O₂
- ◇ 3 times Air as O₂ yields ~ 40% O₂
- ◇ 3 times O₂ as Air yields ~ 80% O₂



Hook-Up Supplies Included:

- 000A0489- DISS O₂ Supply Hose, FxF 4ft Grn
- 000A2420B - Airway Press. Sampling Tee Assy
- 000A5267 - DISS O₂ Supply Hose, FxF 18 in. Grn
- 000A6546 - DISS Air x O₂ Supply Hose, F x F 4ft. Ylw
- 150A1691 - O₂, DISS HT F Inlet x DISS Demand Outlets Y-Block (Grn)
- 150A6016 - O₂, DISS F x 1/4" MNPT x 3/4" L, 90° Elbow
- 240A1455 - Humidifier Bubbler, Refillable
- 240A1457 - Exhalation Valve
- 200A2682 - Wye, Parallel 22mmOD x 22mmOD x 22mmOD/15mmID
- 202A2895 - Tube, gray BS w/cuffs 22mm x 36"
- 202A2896 - Tube, gray BS w/cuffs 22mm x 48" (2)